

## Product datasheet for CP565661

### Tissue Protein Lysate, Kidney

#### Product data:

Product Type:	Tissue Protein Lysates
Disease State:	Other Disease
Diagnosis Category:	Polycystic Kidney Disease
Tissue:	Kidney
Frozen Sample ID:	FR0003BA6F
Case ID:	CI0000014335
Prep:	1P
Vial ID:	PL000027D2
Tissue of Origin:	Kidney
Site of Finding:	Kidney
Appearance:	Lesional
Sample Diagnosis (from Pathology Verification):	Polycystic kidney disease
Normal:	0%
Lesional:	100%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	No normal architecture in lesional component
CASE Diagnosis (from medical center path report):	Polycystic kidney disease
Age:	64



[View online »](#)

<b>Gender:</b>	Male
<b>Race:</b>	White or Caucasian
<b>Tumor Grade:</b>	Not applicable
<b>Minimum Stage Grouping:</b>	Not applicable
<b>T (Extent of Primary Tumor):</b>	Not applicable
<b>N (Lymph node metastasis):</b>	Not applicable
<b>M (Distant metastasis):</b>	Not applicable
<b>Storage:</b>	Store the tissue lysate at -80°C.
<b>Stability:</b>	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.