

## Product datasheet for CP565667

### Tissue Protein Lysate, Lung

#### Product data:

Product Type:	Tissue Protein Lysates
Disease State:	Cancer
Diagnosis Category:	Lung Cancer
Tissue:	Lung
Frozen Sample ID:	FR0002FFCE
Case ID:	CI0000015032
Prep:	2P
Vial ID:	PL000021B2
Tissue of Origin:	Lung
Site of Finding:	Lung
Appearance:	Tumor
Sample Diagnosis (from Pathology Verification):	Adenocarcinoma of lung, acinar, papillary
Normal:	5%
Lesional:	0%
Tumor:	80%
Tumor Hypo/Acellular Stroma:	15%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	Inflammation: Mild Lymphoplasmacytic infiltrate
CASE Diagnosis (from medical center path report):	Adenocarcinoma of lung, acinar, papillary
Age:	77



[View online »](#)

<b>Gender:</b>	Female
<b>Race:</b>	White or Caucasian
<b>Tumor Grade:</b>	AJCC G2: Moderately differentiated
<b>Minimum Stage Grouping:</b>	IIB
<b>T (Extent of Primary Tumor):</b>	T2
<b>N (Lymph node metastasis):</b>	N1
<b>M (Distant metastasis):</b>	MX
<b>Diagnostic Test Results:</b>	ERBB1 ~ EGFR ~ Epidermal growth factor receptor ! by IHC,Positive, moderate. HER2 ~ ERBB2 ~ Neu ! by IHC,Negative
<b>Storage:</b>	Store the tissue lysate at -80°C.
<b>Stability:</b>	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.