

## Product datasheet for CP565683

### Tissue Protein Lysate, Spleen

#### Product data:

Product Type:	Tissue Protein Lysates
Disease State:	Normal
Diagnosis Category:	Normal
Tissue:	Spleen
Frozen Sample ID:	FR71B9C9ED
Case ID:	CI0000016479
Prep:	1P
Vial ID:	PL000022E3
Tissue of Origin:	Spleen
Site of Finding:	Spleen
Appearance:	Normal
Sample Diagnosis (from Pathology Verification):	Within normal limits
Normal:	100%
Lesional:	0%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	Non Tumor Structures: 65% Red pulp, 35% White pulp
CASE Diagnosis (from medical center path report):	Adenocarcinoma of ovary, serous
Age:	73



[View online »](#)

<b>Gender:</b>	Female
<b>Race:</b>	White or Caucasian
<b>Tumor Grade:</b>	Not applicable
<b>Minimum Stage Grouping:</b>	Not applicable
<b>T (Extent of Primary Tumor):</b>	Not applicable
<b>N (Lymph node metastasis):</b>	Not applicable
<b>M (Distant metastasis):</b>	Not applicable
<b>Diagnostic Test Results:</b>	Pancreatic Polypeptide ~ PP,Positive. Estrogen receptor ~ ER ! by IHC,Positive
<b>Storage:</b>	Store the tissue lysate at -80°C.
<b>Stability:</b>	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.