

## Product datasheet for CP565750

### Tissue Protein Lysate, Kidney

#### Product data:

|   |  |
|---|--|
| Product Type:                                     | Tissue Protein Lysates   |
| Disease State:                                    | Other Disease  |
| Diagnosis Category:                               | Polycystic Kidney Disease  |
| Tissue:   | Kidney   |
| Frozen Sample ID:                                 | FR000040E3   |
| Case ID:  | CU0000001295   |
| Prep:   | 3P   |
| Vial ID:  | PL000027D4   |
| Tissue of Origin:                                 | Kidney   |
| Site of Finding:                                  | Kidney   |
| Appearance:                                       | Lesional   |
| Sample Diagnosis (from Pathology Verification):   | Polycystic kidney disease  |
| Normal:   | 0%   |
| Lesional:   | 100%   |
| Tumor:  | 0%   |
| Tumor Hypo/Acellular Stroma:                      | 0%   |
| Tumor Hypercellular Stroma:                       | 0%   |
| Necrosis:   | 0%   |
| Pathology Verification notes from H&E review:     | 45% Cysts with calcifications, 55% Dense interstitial fibrosis present with chronic inflammation |
| CASE Diagnosis (from medical center path report): | Polycystic kidney disease  |
| Age:  | 46   |



[View online »](#)

|                                     |  |
|-------------------------------------|--|
| <b>Gender:</b>                      | Female   |
| <b>Tumor Grade:</b>                 | Not applicable   |
| <b>Minimum Stage Grouping:</b>      | Not applicable   |
| <b>T (Extent of Primary Tumor):</b> | Not applicable   |
| <b>N (Lymph node metastasis):</b>   | Not applicable   |
| <b>M (Distant metastasis):</b>      | Not applicable   |
| <b>Storage:</b>                     | Store the tissue lysate at -80°C.  |
| <b>Stability:</b>                   | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |