

## Product datasheet for CP565823

### Tissue Protein Lysate, Kidney

#### Product data:

|                                                   |                                              |
|---------------------------------------------------|----------------------------------------------|
| Product Type:                                     | Tissue Protein Lysates                       |
| Disease State:                                    | Other Disease                                |
| Diagnosis Category:                               | Polycystic Kidney Disease                    |
| Tissue:                                           | Kidney                                       |
| Frozen Sample ID:                                 | FR00032BAC                                   |
| Case ID:                                          | CU0000011183                                 |
| Prep:                                             | 1P                                           |
| Vial ID:                                          | PL000027D9                                   |
| Tissue of Origin:                                 | Kidney                                       |
| Site of Finding:                                  | Kidney                                       |
| Appearance:                                       | Lesional                                     |
| Sample Diagnosis (from Pathology Verification):   | Polycystic kidney disease                    |
| Normal:                                           | 0%                                           |
| Lesional:                                         | 100%                                         |
| Tumor:                                            | 0%                                           |
| Tumor Hypo/Acellular Stroma:                      | 0%                                           |
| Tumor Hypercellular Stroma:                       | 0%                                           |
| Necrosis:                                         | 0%                                           |
| Pathology Verification notes from H&E review:     | No normal architecture in lesional component |
| CASE Diagnosis (from medical center path report): | Polycystic kidney disease                    |
| Age:                                              | 36                                           |


[View online »](#)

|                              |                                                                                                                  |
|------------------------------|------------------------------------------------------------------------------------------------------------------|
| Gender:                      | Male                                                                                                             |
| Tumor Grade:                 | Not applicable                                                                                                   |
| Minimum Stage Grouping:      | Not applicable                                                                                                   |
| T (Extent of Primary Tumor): | Not applicable                                                                                                   |
| N (Lymph node metastasis):   | Not applicable                                                                                                   |
| M (Distant metastasis):      | Not applicable                                                                                                   |
| Storage:                     | Store the tissue lysate at -80°C.                                                                                |
| Stability:                   | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |