

Product datasheet for CP565843

Tissue Protein Lysate, Lung

Product data:

Product Type:	Tissue Protein Lysates
Disease State:	Normal
Diagnosis Category:	Normal
Tissue:	Lung
Frozen Sample ID:	FR00039C3F
Case ID:	CU0000014403
Prep:	1P
Vial ID:	PL00002327
Tissue of Origin:	Lung
Site of Finding:	Lung
Appearance:	Normal
Sample Diagnosis (from Pathology Verification):	Within normal limits
Normal:	100%
Lesional:	0%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	Non Tumor Structures: 40% Alveoli, 20% Fibrovascular septa, 40% Intrapulmonary lymph node
CASE Diagnosis (from medical center path report):	Adenocarcinoma of lung
Age:	46



[View online »](#)

Gender:	Male
Tumor Grade:	Not applicable
Minimum Stage Grouping:	Not applicable
T (Extent of Primary Tumor):	Not applicable
N (Lymph node metastasis):	Not applicable
M (Distant metastasis):	Not applicable
Storage:	Store the tissue lysate at -80°C.
Stability:	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.