

## Product datasheet for CR559206

### Tissue Total RNA, Liver

#### Product data:

|   |                         |
|---|-------------------------|
| Product Type:                                   | Total RNA               |
| Disease State:                                  | Normal                  |
| Diagnosis Category:                             | Normal                  |
| Tissue:   | Liver                   |
| Frozen Sample ID:                               | FR65EE1885              |
| Case ID:  | CU0000006243            |
| Prep:   | 6                       |
| Vial ID:  | RN00001472              |
| Tissue of Origin:                               | Liver                   |
| Site of Finding:                                | Liver                   |
| Appearance:                                     | Normal                  |
| Sample Diagnosis (from Pathology Verification): | Within normal limits    |
| Normal:   | 100%                    |
| Lesional:                                       | 0%                      |
| Tumor:  | 0%                      |
| Tumor Hypo/Acellular Stroma:                    | 0%                      |
| Tumor Hypercellular Stroma:                     | 0%                      |
| Necrosis:                                       | 0%                      |
| Pathology Verification notes from H&E review:   | 90% lobules, 10% triads |
| Bioanalyzer Ratio (28S/18S):                    | 1.1                     |



[View online »](#)

|  |  |
|--|--|
| <b>CASE Diagnosis (from medical center path report):</b> | Granuloma of liver   |
| <b>Age:</b>  | 52   |
| <b>Gender:</b>   | Male   |
| <b>Tumor Grade:</b>                                      | Not applicable   |
| <b>Minimum Stage Grouping:</b>                           | Not applicable   |
| <b>T (Extent of Primary Tumor):</b>                      | Not applicable   |
| <b>N (Lymph node metastasis):</b>                        | Not applicable   |
| <b>M (Distant metastasis):</b>                           | Not applicable   |
| <b>Storage:</b>  | Store tissue total RNA at -80°C.   |
| <b>Stability:</b>  | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |