

## Product datasheet for CR560421

### Tissue Total RNA, Esophagus

#### Product data:

|   |   |
|---|---|
| Product Type:                                   | Total RNA   |
| Disease State:                                  | Normal  |
| Diagnosis Category:                             | Normal  |
| Tissue:   | Esophagus   |
| Frozen Sample ID:                               | FR5B341F1A  |
| Case ID:  | CI0000012663  |
| Prep:   | 2   |
| Vial ID:  | RN00001F35  |
| Tissue of Origin:                               | Esophagus   |
| Site of Finding:                                | Esophagus   |
| Appearance:                                     | Normal  |
| Sample Diagnosis (from Pathology Verification): | Within normal limits  |
| Normal:   | 100%  |
| Lesional:                                       | 0%  |
| Tumor:  | 0%  |
| Tumor Hypo/Acellular Stroma:                    | 0%  |
| Tumor Hypercellular Stroma:                     | 0%  |
| Necrosis:                                       | 0%  |
| Pathology Verification notes from H&E review:   | Inflammation: Mild Lymphoid infiltrate; Non Tumor Structures: 35% Mucosa, 55% Submucosa, 10% Muscle |
| Bioanalyzer Ratio (28S/18S):                    | 1.55  |



[View online »](#)

|  |  |
|--|--|
| <b>CASE Diagnosis (from medical center path report):</b> | Adenocarcinoma of esophagus  |
| <b>Age:</b>  | 66   |
| <b>Gender:</b>   | Female   |
| <b>Race:</b>   | White or Caucasian   |
| <b>Tumor Grade:</b>                                      | Not applicable   |
| <b>Minimum Stage Grouping:</b>                           | Not applicable   |
| <b>T (Extent of Primary Tumor):</b>                      | Not applicable   |
| <b>N (Lymph node metastasis):</b>                        | Not applicable   |
| <b>M (Distant metastasis):</b>                           | Not applicable   |
| <b>Storage:</b>  | Store tissue total RNA at -80°C.   |
| <b>Stability:</b>  | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |