

Product datasheet for CR562399

Tissue Total RNA, Lymph node

Product data:

| | |
|---|---|
| Product Type: | Total RNA |
| Disease State: | Normal |
| Diagnosis Category: | Normal |
| Tissue: | Lymph node |
| Frozen Sample ID: | FR00000AD3 |
| Case ID: | CU0000000300 |
| Prep: | 3 |
| Vial ID: | RN000000D3 |
| Tissue of Origin: | Lymph node |
| Site of Finding: | Lymph node |
| Appearance: | Normal |
| Sample Diagnosis (from Pathology Verification): | Within normal limits |
| Normal: | 100% |
| Lesional: | 0% |
| Tumor: | 0% |
| Tumor Hypo/Acellular Stroma: | 0% |
| Tumor Hypercellular Stroma: | 0% |
| Necrosis: | 0% |
| Pathology Verification notes from H&E review: | Sample contains well formed normal lymph node follicles |
| Bioanalyzer Ratio (28S/18S): | 1.46 |



[View online »](#)

| | |
|--|--|
| CASE Diagnosis (from medical center path report): | Adenocarcinoma of lung, bronchioloalveolar |
| Age: | 41 |
| Gender: | Female |
| Tumor Grade: | Not applicable |
| Minimum Stage Grouping: | Not applicable |
| T (Extent of Primary Tumor): | Not applicable |
| N (Lymph node metastasis): | Not applicable |
| M (Distant metastasis): | Not applicable |
| Storage: | Store tissue total RNA at -80°C. |
| Stability: | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |