

## Product datasheet for **CS500489**

### Frozen Tissue Sections, Myometrium


#### Product data:

|   |   |
|---|---|
| Product Type:                                     | Frozen Sections                             |
| Disease State:                                    | Other Disease                               |
| Diagnosis Category:                               | Uterine Fibroids                            |
| Tissue:   | Myometrium                                  |
| Frozen Sample ID:                                 | FR00002AE3                                  |
| Case ID:  | CI0000000066                                |
| Tissue of Origin:                                 | Myometrium                                  |
| Site of Finding:                                  | Myometrium                                  |
| Appearance:                                       | T   |
| Sample Diagnosis (from Pathology Verification):   | Leiomyoma of myometrium                     |
| Normal:   | 0%  |
| Lesional:   | 0%  |
| Tumor:  | 100%  |
| Tumor Hypo/Acellular Stroma:                      | 0%  |
| Tumor Hypercellular Stroma:                       | 0%  |
| Necrosis:   | 0%  |
| Pathology Verification notes from H&E review:     | Leiomyoma                                   |
| CASE Diagnosis (from medical center path report): | Adenocarcinoma of endometrium, endometrioid |
| Age:  | 62  |
| Gender:   | Female                                      |
| Race:   | Black or African American                   |




[View online »](#)

|                              |  |
|------------------------------|--|
| Tumor Grade:                 | FIGO G3: Poorly differentiated   |
| Minimum Stage Grouping:      | IIB  |
| T (Extent of Primary Tumor): | T2b  |
| N (Lymph node metastasis):   | N0   |
| M (Distant metastasis):      | MX   |
| Storage:                     | Store frozen tissue sections at -80°C.   |
| Stability:                   | All frozen tissue sections are stable for 1 year from date of purchase if stored at -80°C without repeated freeze/thaws. |

**Product images:**

4x Image



20x Image