

Product datasheet for CS504067

Frozen Tissue Sections, Kidney

Product data:

Product Type:	Frozen Sections
Disease State:	Other Disease
Diagnosis Category:	Nephritis
Tissue:	Kidney
Frozen Sample ID:	FR000082B7
Case ID:	CI0000000504
Tissue of Origin:	Kidney
Site of Finding:	Kidney
Appearance:	L
Sample Diagnosis (from Pathology Verification):	Pyelonephritis, chronic
Normal:	0%
Lesional:	100%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	10% cortex, 90% medulla
CASE Diagnosis (from medical center path report):	Pyelonephritis, chronic
Age:	32
Gender:	Female
Race:	Not reported




[View online »](#)

Diagnostic Test Results: PAS ~ Periodic acid-Schiff stain, Reported as Performed or Reviewed; Von Kossa, Reported as Performed or Reviewed; Hale colloidal iron stain, Reported as Performed or Reviewed


Storage: Store frozen tissue sections at -80°C.

Stability: All frozen tissue sections are stable for 1 year from date of purchase if stored at -80°C without repeated freeze/thaws.

Product images:



4x Image



20x Image