

## Product datasheet for **CS505683**

### Frozen Tissue Sections, Colon

#### Product data:

|   |                                     |
|---|-------------------------------------|
| Product Type:                                     | Frozen Sections                     |
| Disease State:                                    | Cancer                              |
| Diagnosis Category:                               | Colorectal Cancer                   |
| Tissue:   | Colon                               |
| Frozen Sample ID:                                 | FR00019C75                          |
| Case ID:  | CI0000005399                        |
| Tissue of Origin:                                 | Colon                               |
| Site of Finding:                                  | Liver                               |
| Appearance:                                       | T                                   |
| Sample Diagnosis (from Pathology Verification):   | Adenocarcinoma of colon, metastatic |
| Normal:   | 10%                                 |
| Lesional:   | 0%                                  |
| Tumor:  | 55%                                 |
| Tumor Hypo/Acellular Stroma:                      | 28%                                 |
| Tumor Hypercellular Stroma:                       | 5%                                  |
| Necrosis:   | 2%                                  |
| CASE Diagnosis (from medical center path report): | Adenocarcinoma of colon, metastatic |
| Age:  | 67                                  |
| Gender:   | Male                                |
| Race:   | Not reported                        |
| Tumor Grade:                                      | AJCC G2: Moderately differentiated  |



[View online »](#)

Minimum Stage Grouping: IV

T (Extent of Primary Tumor): TX


N (Lymph node metastasis): NX

M (Distant metastasis): M1


Storage: Store frozen tissue sections at -80°C.

Stability: All frozen tissue sections are stable for 1 year from date of purchase if stored at -80°C without repeated freeze/thaws.

### Product images:



4x Image



20x Image