

Product datasheet for **CS523977**

Frozen Tissue Sections, Tendon

Product data:

Product Type:	Frozen Sections
Disease State:	Normal
Diagnosis Category:	Normal
Tissue:	Tendon
Frozen Sample ID:	FR0002BE93
Case ID:	CI0000009267
Tissue of Origin:	Tendon
Site of Finding:	Tendon
Appearance:	R
Sample Diagnosis (from Pathology Verification):	Within normal limits
Normal:	100%
Lesional:	0%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	100% tendon
CASE Diagnosis (from medical center path report):	Osteoarthritis
Age:	70
Gender:	Female
Race:	Black or African American




[View online »](#)


Storage: Store frozen tissue sections at -80°C .

Stability: All frozen tissue sections are stable for 1 year from date of purchase if stored at -80°C without repeated freeze/thaws.

Product images:



4x Image



20x Image