

## Product datasheet for **CS535825**

### Frozen Tissue Sections, Ovary: left

#### Product data:

|   |                                      |
|---|--------------------------------------|
| Product Type:                                     | Frozen Sections                      |
| Disease State:                                    | Other Disease                        |
| Diagnosis Category:                               | Borderline Ovarian Tumor             |
| Tissue:   | Ovary                                |
| Frozen Sample ID:                                 | FR00039799                           |
| Case ID:  | CI0000014167                         |
| Tissue of Origin:                                 | Ovary: left                          |
| Site of Finding:                                  | Ovary: left                          |
| Appearance:                                       | T                                    |
| Sample Diagnosis (from Pathology Verification):   | Tumor of ovary, mucinous, borderline |
| Normal:   | 0%                                   |
| Lesional:   | 0%                                   |
| Tumor:  | 40%                                  |
| Tumor Hypo/Acellular Stroma:                      | 40%                                  |
| Tumor Hypercellular Stroma:                       | 20%                                  |
| Necrosis:   | 0%                                   |
| CASE Diagnosis (from medical center path report): | Tumor of ovary, mucinous, borderline |
| Age:  | 63                                   |
| Gender:   | Female                               |
| Race:   | White or Caucasian                   |
| Tumor Grade:                                      | AJCC GB: Borderline malignancy       |



[View online »](#)

Minimum Stage Grouping: IA

T (Extent of Primary Tumor): T1a


N (Lymph node metastasis): NX

M (Distant metastasis): MX


Storage: Store frozen tissue sections at -80°C.

Stability: All frozen tissue sections are stable for 1 year from date of purchase if stored at -80°C without repeated freeze/thaws.

**Product images:**



4x Image



20x Image