

Product datasheet for CS536535

Frozen Tissue Sections, Kidney

Product data:

Product Type:	Frozen Sections
Disease State:	Other Disease
Diagnosis Category:	Polycystic Kidney Disease
Tissue:	Kidney
Frozen Sample ID:	FR0003BA73
Case ID:	CI0000014335
Tissue of Origin:	Kidney
Site of Finding:	Kidney
Appearance:	L
Sample Diagnosis (from Pathology Verification):	Polycystic kidney disease
Normal:	0%
Lesional:	100%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	No normal architecture in lesional component
CASE Diagnosis (from medical center path report):	Polycystic kidney disease
Age:	64
Gender:	Male
Race:	White or Caucasian




[View online »](#)


Storage: Store frozen tissue sections at -80°C.

Stability: All frozen tissue sections are stable for 1 year from date of purchase if stored at -80°C without repeated freeze/thaws.

Product images:



4x Image



20x Image