

Product datasheet for **CS543635**

Frozen Tissue Sections, Tendon

Product data:

Product Type:	Frozen Sections
Disease State:	Other Disease
Diagnosis Category:	Arthritis
Tissue:	Tendon
Frozen Sample ID:	FR00040E47
Case ID:	CI0000017873
Tissue of Origin:	Tendon
Site of Finding:	Tendon
Appearance:	L
Sample Diagnosis (from Pathology Verification):	Arthritis, rheumatoid
Normal:	0%
Lesional:	100%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	Non Tumor Structures: 100% Tendon
CASE Diagnosis (from medical center path report):	Arthritis, rheumatoid
Age:	61
Gender:	Female
Race:	White or Caucasian




[View online »](#)


Storage: Store frozen tissue sections at -80°C.

Stability: All frozen tissue sections are stable for 1 year from date of purchase if stored at -80°C without repeated freeze/thaws.

Product images:



4x Image



20x Image