

## Product datasheet for **CS815506**

### PD1 Positive IHC Control Slides


#### Product data:

Product Type: IHC Control Slides (Phase II)


Name: PD1 Positive

Description: PD1 positive tissue 4um section mounted on unstained slide. Available in sets of 2 or 5 slides.

#### Product images:



Hematoxylin and eosin stain of normal tonsil.



Immunohistochemical staining of FFPE tonsil - normal adjacent; using heat-induced epitope retrieval HIER ACCEL buffer pH8.7 at 120C for 3min, mouse monoclonal antibody CAT# - CE00059 anti-PD1 clone UMAB198 was used at 1:200. Strong staining seen in germinal and non-germinal center of the tonsil.



[View online »](#)