

## Product datasheet for **CV900015L**

### AAV5 Scrambled shRNA Control Particle


#### Product data:

Product Type:	AAV Control Particles
Description:	AAV5 control particles, shRNA scramble, expressing GFP and puro, >10 <sup>13</sup> GC/mL, 1 x 100 uL
Reporter:	GFP
Promoter:	CMV
Serotype:	AAV-5
Purification Method:	Iodixanol
Storage Buffer:	PBS with 0.001% Pluronic F68
Validation:	To validate the quality of our AAV control particles, we have developed the following procedure: <ol style="list-style-type: none"><li>1. AAV genome titer is determined by qPCR.</li><li>2. AAV purity is determined by Silver staining.</li><li>3. AAV transduction efficiency is analyzed by in vitro transduction in HEK293T.</li></ol>
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.



[View online »](#)

**Product images:**



In-vitro transduction efficiency of HEK293T cell line with AAV-5 Scrambled shRNA Control Particles.