

Product datasheet for CV900017S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

AAV8 Control Particle

Product data:

Product Type: AAV Control Particles

Description: AAV8 control particles, shRNA scramble, expressing GFP and puro, >10^13 GC/mL, 50ul

Reporter: GFP Promoter: CMV

Serotype: AAV-8

Purification Method: lodixanol

Storage Buffer: PBS with 0.001% Pluronic F68

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at

-20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, vectors should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.