

## Product datasheet for **CV900019S**

### AAV5 Control Particle

#### Product data:

|                      |   |
|----------------------|---|
| Product Type:        | AAV Control Particles   |
| Description:         | AAV5 with CMV promoter-driven expression of Myc-DDK, >10 <sup>13</sup> GC/mL, 50ul  |
| Reporter:            | Myc-DDK   |
| Promoter:            | CMV   |
| Serotype:            | AAV-5   |
| Purification Method: | Iodixanol   |
| Storage Buffer:      | PBS with 0.001% Pluronic F68  |
| Stability:           | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, vectors should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |



[View online »](#)