

## Product datasheet for **CV900021L**

### AAV8-CMV-Myc-DDK control particles

#### Product data:

Product Type:	AAV Control Particles
Description:	AAV8 with CMV promoter-driven expression of Myc-DDK, $>10^{13}$ GC/mL, 1 x 100 uL
Reporter:	Myc-DDK
Promoter:	CMV
Serotype:	AAV-8
Purification Method:	Iodixanol
Storage Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, vectors should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.



[View online »](#)