


## Product datasheet for **CV900023S**

### AAV1 Scrambled shRNA Control Particle

#### Product data:

Product Type:	AAV Control Particles
Description:	AAV1 control particles, shRNA scramble, expressing GFP and puro, >10 <sup>13</sup> GC/mL, 50 ul
Reporter:	GFP
Promoter:	CMV
Serotype:	AAV-1
Purification Method:	Iodixanol
Storage Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, vectors should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

#### Product images:



AAV1 Scrambled shRNA Control 10x



[View online »](#)