

Product datasheet for **DA3551X**

PDGFA (PDGF-AA) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	PDGFA (PDGF-AA) human recombinant protein, 20 µg
Species:	Human
Expression Host:	E. coli
Predicted MW:	28.5 kDa
Purity:	>95% pure by SDS-PAGE and visualised by silver stain.
Buffer:	Presentation State: Purified State: Lyophilized purified fraction. Buffer System: 50 mM Acetic Acid
Bioactivity:	Biological: The ED50 as determined by the dose-dependent stimulation of thymidine uptake by BALB/c 3T3 cells is < 1 ng/ml. Specific: > 1 x 10e6 units/mg.
Endotoxin:	< 0.1 ng per µg of PDGF-AA.
Preparation:	Lyophilized purified fraction.
Protein Description:	Recombinant Human PDGF-AA is a 28.5 kDa disulfide-linked homodimer of two A chains (250 total amino acids). Range: 0.2-10.0 ng/ml
Note:	Always centrifuge product before opening vial!
Storage:	Store lyophilized human PDGF-AA at -20°C. Reconstituted PDGF-AA should be stored in working aliquots at -20°C to -70°C. Avoid repeated freeze-thaw cycles!
Stability:	Shelf life: One year from despatch.
RefSeq:	NP_002598
Locus ID:	5154
UniProt ID:	P04085
Cytogenetics:	7p22.3
Synonyms:	PDGF-A; PDGF1



[View online »](#)

Summary:	This gene encodes a member of the protein family comprised of both platelet-derived growth factors (PDGF) and vascular endothelial growth factors (VEGF). The encoded preproprotein is proteolytically processed to generate platelet-derived growth factor subunit A, which can homodimerize, or alternatively, heterodimerize with the related platelet-derived growth factor subunit B. These proteins bind and activate PDGF receptor tyrosine kinases, which play a role in a wide range of developmental processes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2015]
Protein Families:	Druggable Genome
Protein Pathways:	Cytokine-cytokine receptor interaction, Focal adhesion, Gap junction, Glioma, MAPK signaling pathway, Melanoma, Pathways in cancer, Prostate cancer, Regulation of actin cytoskeleton