


Product datasheet for **DDX0200A546-100**

DC SIGN (CD209) (Extracell. Dom.) Mouse Monoclonal Antibody [Clone ID: 118A8.05]

Product data:

Product Type:	Primary Antibodies
Clone Name:	118A8.05
Applications:	IF
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	Label: Alexa Fluor 546
Concentration:	0.5 mg/ml
Conjugation:	Alexa Fluor 546
Gene Name:	CD209 molecule
Database Link:	Entrez Gene 30835 Human
Synonyms:	DCSIGN1, DCSIGN, DC-SIGN1, CLEC4L, Dendritic Cell Marker

Product images:



IHC staining of human skin frozen section with clone 118A8 (DX0200)



[View online »](#)