

Product datasheet for **DDX0451P-100**

Dystonin (DST) Human Monoclonal Antibody [Clone ID: H10D-Hy-8B]

Product data:

Product Type: Primary Antibodies
Clone Name: H10D-Hy-8B
Applications: IF, IHC, WB
Recommend Dilution: **WB, IF, IHC frozen section.**

Usage recommendation:

*This monoclonal antibody may be used between 5-20 µg/ml.
*Optimal dilution should be determined by each laboratory for each application.
*Coupled antibody: to maintain RT before using.

Reactivity: Human
Host: Human
Isotype: IgG4
Clonality: Monoclonal
Immunogen: EBV cell clones derived from PBMC of patients with the typical clinical, histologic, immunopathologic criteria of Bullous Pemphigoid
Formulation: **Purified:** 100 µg in 200µl / 50 µg in 100 µl Tris-NaCl pH 8.
Coupled: 100 µg in 200 µl / 50 µg in 100 µl PBS 50% glycerol.
State: Purified
Concentration: 0.5 mg/ml
Purification: QMA Hyper D ion exchange chromatography
Gene Name: dystonin
Database Link: [Entrez Gene 667 Human](#)
Background: Bullous pemphigoid is an acquired blistering skin disease associated with production of IgG autoantibodies to the BP antigen. Antibodies-producing cells were established by EBV infection of CD40-activated PBMCs. These antibodies recognize different epitopes of BPAG1. (Peyron E. et al., J Immunol, 1994, 153, 1333)
Synonyms: Hemidesmosomal plaque protein, Dystonia musculorum protein, BPA, Bullous pemphigoid antigen, DST, BPAG1, DMH, DT, KIAA0728



[View online »](#)