

Product datasheet for **DDX0463A488-100**

SEMA6A Mouse Monoclonal Antibody [Clone ID: 118E7]

Product data:

Product Type:	Primary Antibodies
Clone Name:	118E7
Applications:	FC
Recommend Dilution:	<u>DDX0463P-50 / DDX0463P-100 Purified:</u> Surface Flow Cytometry. <u>DDX0463A488-50 / DDX0463A488-100 Alexa-fluor®488:</u> Surface Flow Cytometry. <u>DDX0463A546-50 / DDX0463A546-100 Alexa-fluor®546:</u> IF. <u>DDX0463A647-50 / DDX0463A647-100 Alexa- fluor®647:</u> (on request) IHC, Flow cytometry.

Usage recommendation:

- *This monoclonal antibody may be used between 5-20 µg/ml.
- *Optimal dilution should be determined by each laboratory for each application.
- *Coupled antibody: to maintain RT before using.

Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	human recombinant SEMA-6A
Specificity:	Human SEMA-6A (epitope in extracellular domain).
Formulation:	<u>Purified:</u> 100 µg in 200µl / 50 µg in 100 µl Tris-NaCl pH 8. <u>Coupled:</u> 100 µg in 200µl / 50 µg in 100 µl PBS 50% glycerol. Label: Alexa Fluor 488
Concentration:	0.5 mg/ml
Purification:	QMA Hyper D ion exchange chromatography
Conjugation:	Alexa Fluor 488
Gene Name:	semaphorin 6A
Database Link:	Entrez Gene 57556 Human



[View online »](#)

Background:

Semaphorins represent a large family composed of 19 secreted, GPI-linked or transmembrane members divided into 5 classes (3 to 7) and characterized by the presence of a 500 amino acid Sema domain in the Nterminus. Initially described for their role in axon guidance and apoptosis during neuronal development, some semaphorins like the class 4 regulate immune functions and are involved in organogenesis and angiogenesis. The transmembrane class 6 semaphorins include 4 members (SEMA-6A to 6D). SEMA-6A features an intracellular binding site for regulatory Evi protein. In the immune system, SEMA6A protein is restricted in vitro and in vivo to an activation status of Langerhans cells. (*Gautier G et al., 2006 Am J Pathol., 168: 453-65*).

Synonyms:

Semaphorin-6A, Semaphorin VIA, Sema VIA, Semaphorin-6A-1, SEMA6A-1, SEMAQ, KIAA1368

Product images:
