

Product datasheet for **DDX0492P-100**

TLR10 Mouse Monoclonal Antibody [Clone ID: 9F4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	9F4
Applications:	IHC
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	State: Purified
Gene Name:	toll like receptor 10
Database Link:	Entrez Gene 81793 Human
Synonyms:	Toll-like receptor 10, UNQ315/PRO358



[View online »](#)