

## Product datasheet for DM001

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## alpha smooth muscle Actin (ACTA2) Mouse Monoclonal Antibody [Clone ID: 1A4/asm-1]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: 1A4/asm-1

Applications: IHC

Recommended Dilution: Immunohistochemistry on Formalin Fixed Paraffin Embedded Sections: Use Actin

antibody at 1/25-1/50 dilution in an ABC method for 30 minutes at room temperature.

**Recommended Positive Control:** Leiomyoma colon.

**Reactivity:** Bovine, Chicken, Human, Mouse, Rabbit, Rat

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: N-terminal decapeptide of alpha-smooth muscle Actin

**Specificity:** This antibody is specific to alpha-smooth isoform of actin.

It reacts with smooth muscle cells of vessels and different parenchymes.

The antibody does not cross-react with beta and gamma-cytoplasmic, alpha-sarcomeric and

alpha-myocardial actin isoforms. *Cellular Localization:* Cytoplasmic.

**Formulation:** State: Purified

State: Liquid purified Ig fraction containing Sodium Azide as preservative

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C.

Stability: Shelf life: one year from despatch.

Gene Name: actin, alpha 2, smooth muscle, aorta

Database Link: Entrez Gene 59 Human

P62736






Background:

Smooth muscle actin is expressed in normal smooth muscle and in myoepithelial cells, and also in tumours derived from these cells. Actins are involved in maintaining the structure and integrity of a cell, as well as cell motility. They are highly conserved proteins and make up the majority of the contractile apparatus. To date, six different actin isoforms have been described.


Synonyms:

ACTSA, ACTVS, Alpha-actin-2

## **Product images:**



Formalin fixed paraffin embedded human leiomyoma stained with Smooth Muscle Actin antibody DM001.



Formalin fixed paraffin embedded human colon stained with Smooth Muscle Actin antibody DM001.