

## Product datasheet for **DM052-05**


### Cytokeratin 5+6+8+17+19 Mouse Monoclonal Antibody [Clone ID: MNF116]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	MNF116
Applications:	IHC
Recommended Dilution:	<b>Immunohistochemistry on Paraffin Embedded Sections:</b> Use a dilution of 1/25-1/50 in an ABC method. Enzymatic predigestion of paraffin embedded tissues is recommended prior to the immunostaining. <b>Incubation Time:</b> 30-60 minutes at RT. <b>Positive Control:</b> Skin.
Reactivity:	Bovine, Equine, Guinea Pig, Human, Mouse, Porcine, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Crude extract of splenic cells from a nude mouse engrafted with MCF-7 cells.
Specificity:	This antibody recognizes cytokeratin polypeptides ranging from 40 through 58kD, corresponding to Cytokeratins 5, 6, 8, 17 and 19. It shows a broad pattern of reactivity with Human epithelial tissues from simple glandular epithelia to stratified squamous epithelia, like epidermis, mammary gland ducts and tracheal epithelium. <b>Cellular Localization:</b> Cytoplasmic.
Formulation:	State: Supernatant State: Liquid Tissue Culture Supernatant Preservative: Sodium Azide
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Synonyms:	Cytokeratin Pan-reactive



[View online »](#)

**Product images:**

Formalin fixed paraffin embedded human tonsil stained with Cytokeratin Antibody (Clone MNF116).