

Product datasheet for **DM074**

Human IgM Mouse Monoclonal Antibody [Clone ID: R1/69]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM Mouse Monoclonal Antibody [Clone ID: R1/69]
Clone Name:	R1/69
Applications:	IHC
Recommended Dilution:	Immunohistochemistry on Formalin-Fixed Paraffin Embedded Sections: Use a dilution of 1/25-1/50 in an ABC method (30 minutes at room temperature). Enzymatic predigestion with proteolytic enzymes is recommended for paraffin embedded tissue sections. <i>Recommended Positive Control:</i> Tonsil.
Reactivity:	Human
Host:	Mouse
Immunogen:	Human lymphoid cells.
Isotype:	IgG1
Formulation:	State: Supernatant State: Tissue Culture Supernatant containing Sodium Azide as preservative.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)