

Product datasheet for **DM1032**

Prostate Mucin Antigen / PMA Mouse Monoclonal Antibody [Clone ID: YPMA-2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	YPMA-2
Applications:	IF, IHC, WB
Recommended Dilution:	Immunoperoxidase staining. Immunofluorescent staining. Western Blot.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG3
Clonality:	Monoclonal
Immunogen:	Membrane protein extracted from human prostate mucin tissues
Specificity:	This antibody reacts to prostate mucin antigen. Does not react with fetal or benign prostate specimens, non-prostate carcinomas and normal human tissues.
Formulation:	0.01M PBS, pH 7.2 State: Purified State: Lyophilized purified IgG
Reconstitution Method:	Restore with double distilled water is recommended and to adjust the final concentration to 1.00 mg/ml
Purification:	Protein G affinity purified
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)