

Product datasheet for DM2036

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

MCP3 (CCL7) Mouse Monoclonal Antibody [Clone ID: 343]

Product data:

Product Type: Primary Antibodies

Clone Name: 343

Applications: ELISA, WB Recommended Dilution: ELISA.

Western Blot: An antibody concentration of 1-4 µg/mL will allow visualization of 0.1 µg/lane of

human MCP-3.

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Recombinant human MCP-3 (M.W. 8.5 kDa containing 76 amino acids residues)

Specificity: Reactive with natural and recombinant human MCP-3.

Does not show any cross reaction with other Human Cytokines or growth factors tested such

as IL1 beta, IL-8, MCAF, SAA and EGF

Formulation: 0.01M PBS, pH 7.2 without preservatives.

State: Purified

State: Lyophilized purified IgG fraction.

Reconstitution Method: Restore with double distillated water to adjust the final concentration to 1.00 mg/ml

Purification: Affinity Chromatography on Protein G.

Conjugation: Unconjugated

Storage: Store the antibody at -20°C.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: C-C motif chemokine ligand 7

Database Link: Entrez Gene 6354 Human

P80098



MCP3 (CCL7) Mouse Monoclonal Antibody [Clone ID: 343] - DM2036

Background: MCP3/CCL7 is a secreted chemokine which attracts macrophages during inflammation and

metastasis. It is a member of the C-C subfamily of chemokines which are characterized by

having two adjacent cysteine residues. The protein is an in vivo substrate of matrix metalloproteinase 2, an enzyme which degrades components of the extracellular matrix.

Synonyms: Small-inducible cytokine A7, CCL-7, MCP-3, C-C motif chemokine 7, SCYA6, SCYA7, NC28