

Product datasheet for **DM3526**

PLGF (PGF) Mouse Monoclonal Antibody [Clone ID: 331/H12]

Product data:

Product Type:	Primary Antibodies
Clone Name:	331/H12
Applications:	WB
Recommended Dilution:	Western Blot: 1-2 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant human PIGF-2 (produced in insect cells)
Specificity:	This antibody reacts to PIGF.
Formulation:	PBS, pH 7.4 State: Purified State: Lyophilized purified Ig
Reconstitution Method:	Sterile water to a concentration of 1.0 mg/ml
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch
Gene Name:	placental growth factor
Database Link:	Entrez Gene 5228 Human P49763
Synonyms:	PGFL, PLGF, PIGF



[View online »](#)