

Product datasheet for **DM3540P**

Dhh (C-term) Rat Monoclonal Antibody [Clone ID: 7C8]

Product data:

Product Type:	Primary Antibodies
Clone Name:	7C8
Applications:	WB
Recommended Dilution:	Western Blot: 1/500 - 1/1000.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2
Clonality:	Monoclonal
Immunogen:	Recombinant dhh C-terminal fragment
Specificity:	This antibody detects dhh N-terminal polypeptide, but not C-terminal.
Formulation:	State: Purified State: Lyophilized purified Ig
Reconstitution Method:	Restore with 200µl sterile PBS and the final concentration is 500µg/ml.
Purification:	Protein A/G affinity chromatography
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -70°C. Following reconstitution store the antibody (in aliquots) at -20°C for 6 month. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	desert hedgehog
Database Link:	Entrez Gene 13363 Mouse Q61488
Synonyms:	HHG-3



[View online »](#)