

Product datasheet for **DM3572P**

IL15ra (Extracell. Dom.) Rat Monoclonal Antibody [Clone ID: 6C38]

Product data:

Product Type:	Primary Antibodies
Clone Name:	6C38
Applications:	WB
Recommended Dilution:	Western Blot: 1/500 - 1/1000.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2
Clonality:	Monoclonal
Immunogen:	Recombinant protein of IL-15Ralpha extracellular domain.
Specificity:	This antibody detects IL-15Ralpha. No cross activity to other IGFbps.
Formulation:	State: Purified State: Lyophilized purified Ig
Reconstitution Method:	Restore with 200µl sterile PBS and the final concentration is 500µg/ml.
Purification:	Protein A/G affinity chromatography
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -70°C. Following reconstitution store the antibody (in aliquots) at -20°C for 6 month. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 15 receptor, alpha chain
Database Link:	Entrez Gene 16169 Mouse Q60819



[View online »](#)

Background:

IL15RA is a cytokine receptor that specifically binds IL15 with high affinity. It shares two subunits with the receptor of IL2, the IL2R beta and IL2R gamma chains. This forms the basis of many overlapping biological activities of IL15 and IL2. The IL2 receptor requires an additional IL2-specific alpha subunit for high affinity IL2 binding. This protein is structurally related to IL2R alpha, but is capable of binding IL15 with high affinity independent of other subunits, which suggests the distinct roles between IL15 and IL2. This receptor is reported to enhance cell proliferation and expression of apoptosis inhibitor BCL2L1/BCL2-XL and BCL2. There are nine different isoforms.

Synonyms:

Interleukin-15 receptor subunit alpha, IL-15R-alpha, IL-15RA