

Product datasheet for **DM3575P**

IL18r1 (Extracell. Dom.) Rat Monoclonal Antibody [Clone ID: 7G31]

Product data:

Product Type:	Primary Antibodies
Clone Name:	7G31
Applications:	WB
Recommended Dilution:	Western Blot: 1/500 - 1/1000.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2
Clonality:	Monoclonal
Immunogen:	Recombinant protein of IL-18R extracellular domain.
Specificity:	This antibody detects IL-18R.
Formulation:	State: Purified State: Lyophilized purified Ig
Reconstitution Method:	Restore with 200µl sterile PBS and the final concentration is 500µg/ml.
Purification:	Protein A/G affinity chromatography
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -70°C. Following reconstitution store the antibody (in aliquots) at -20°C for 6 month. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 18 receptor 1
Database Link:	Entrez Gene 16182 Mouse Q61098
Synonyms:	Interleukin-18 receptor 1, IL-18 receptor 1, CDw218a, IL-1Rrp, IL1RRP



[View online »](#)