


## Product datasheet for DP020

### Human IgA Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA Rabbit Polyclonal Antibody
Applications:	IHC
Recommended Dilution:	Immunohistochemistry on Frozen and Paraffin Sections (ABC method: 1/25-1/100). Formalin-Fixed Paraffin Embedded Tissue Sections requires High temperature antigen unmasking with 10mM citrate buffer, pH 6.0 must be performed prior to the immunostaining. <i>Incubation Time:</i> 30 minutes at room temperature. <i>Positive Control:</i> Tonsil.
Reactivity:	Human
Host:	Rabbit
Immunogen:	IgA isolated from a pool of normal Human sera.
Formulation:	State: Purified State: Liquid purified Ig fraction containing Sodium Azide as preservative.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.

#### Product images:



Formalin fixed paraffin embedded human tonsil stained with IgA antibody. DP020



[View online »](#)