


## Product datasheet for DP023-05

### Human IgG Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG Rabbit Polyclonal Antibody
Applications:	IHC
Recommended Dilution:	<b>Immunohistochemistry on Formalin-Fixed Paraffin Embedded Sections:</b> Use at 1/200-1/500 dilution in an ABC method (30 minutes at room temperature). Proteolytic treatment is required prior to immunostaining of formalin fixed paraffin embedded tissue sections. <i>Recommended Positive Control:</i> Tonsil.
Reactivity:	Human
Host:	Rabbit
Immunogen:	IgG isolated from a pool of normal Human sera.
Formulation:	State: Purified State: Liquid Ig fraction containing Sodium Azide as preservative.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.

#### Product images:



Formalin fixed paraffin embedded human tonsil stained with IgG antibody. DP023



[View online »](#)