

## Product datasheet for **DP024-05**

### Human IgM Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM Rabbit Polyclonal Antibody
Applications:	IHC
Recommended Dilution:	Immunohistochemistry with normal and neoplastic tissues (formalin fixed paraffin-embedded): the antibody may be diluted 1:150-1:250 in an ABC method, incubation time 30 minutes. Positive control: Tonsil. Cellular localization: Cytoplasmic.
Reactivity:	Human
Host:	Rabbit
Immunogen:	IgM isolated from human sera of a patient
Formulation:	State: Purified State: Liquid purified Ig-fraction containing sodium azide as a preservative.
Conjugation:	Unconjugated
Storage:	Store the antibody at 2-8°C.



[View online »](#)