

Product datasheet for DP026-05

Floudet datasileet for DF020-03

Ki67 (MKI67) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC

Recommended Dilution: Immunohistochemsistry on Frozen and Formalin-Fixed Paraffin Embedded Sections:

Use a dilution of 1/25-1/50 in an ABC method (30 minutes at room temperature). Formalin fixed paraffin embedded tissue sections require high temperature antigen

unmasking with 10 mM citrate buffer, pH 6.0 prior to immunostaining.

Recommended Positive Control: Tonsil.

Reactivity: Human

Host: Rabbit

Clonality: Polyclonal

Immunogen: Synthetic peptide from 62 base pair region of the Human Ki-67 antigen.

Specificity: This antibody reacts with a nuclear antigen present in proliferating human cells.

It can be used to evaluate the ki-67 labeling index in various tumors.

Cellular Localization: Nuclear.

Formulation: State: Purified

State: Liquid purified Ig fraction Preservative: < 0.1% Sodium Azide

Conjugation: Unconjugated

Storage: Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: marker of proliferation Ki-67

Database Link: Entrez Gene 4288 Human

P46013



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Background:

Ki-67 antigen is expressed in proliferating cells but not in quiescent cells. Expression of this antigen occurs preferentially during late G1, S, G2, and M phases of the cell cycle, while in cells in G0 phase the antigen cannot be detected. Consequently, Ki-67 antigen expression is used in tumor pathology to detect proliferating cells in neoplastic diseases. In cultured cells, Ki-67 is expressed in the nucleolus of interphase cells. The Ki-67 gene contains 15 exons. The Ki-67 repeat region, within which there is a 22-amino acid Ki-67 motif, is encoded by exon 13. The shorter isoform lacks exon 7. Northern blot analysis re-vealed multiple transcripts ranging from approximately 8.9 to 12.5 kb in proliferating but not quies-cent cells. Immunoblot analysis showed expression of 320- and 359-kD proteins. Sequence analy-sis predicted that the short-lived 2,896- and 3,256-amino acid protein isoforms contain potential nuclear targeting signals, over 200 potential phosphorylation sites, 19 N-myristoylation sites, 3 amidation sites, and numerous PEST sites. Antisense oligonucleotides inhibited cellular prolifera-tion in a dose-dependent manner, suggesting that Ki-67 protein expression may be an absolute requirement for cell proliferation.

Synonyms:

Ki67 antigen, MKI67, Proliferation marker

Product images:



Formalin fixed paraffin embedded human tonsil stained with Ki-67 antibody. DP026