

Product datasheet for **DP3523**

LIF Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Purified
Conjugation:	Unconjugated
Gene Name:	leukemia inhibitory factor
Database Link:	Entrez Gene 3976 Human P15018
Synonyms:	LIF, HILDA, D factor, Melanoma-derived LPL inhibitor, MLPLI, Emfilermin



[View online »](#)