

Product datasheet for **EA101826**

IL2

Product data:

Product Type:	ELISA Kits
Description:	Murine IL-2 ELISpot Pair, 10 x 96-well
Size:	10 x 96 wells
Format:	10 x 96-well
Assay Type:	Enzyme-Linked ImmunoSpot (ELISpot)
Signal:	Colorimetric
Specificity:	Recognizes natural murine IL-2
Cross Reactivity:	No cross reactivity with other murine cytokines
Components:	ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
Background:	cytokine produced by T-cells in response to antigen or mitogen stimulation [RGD, Feb 2006]. Publication Note: This RefSeq record includes a subset of the publications that are available for this gene. Please see the Gene record to access additional publications. ##Evidence-Data-START#Transcript exon combination :: M22899.1 [ECO:0000332] RNAseq introns :: single sample supports all introns SRS369737 [ECO:0000348] ##Evidence-Data-END##
Gene Symbol:	IL2
Gene ID:	116562



[View online »](#)