

## Product datasheet for EA101827

### Ifng

#### Product data:

Product Type:	ELISA Kits
Description:	Murine IFN gamma ELISpot Pair, 10 x 96-well
Size:	10 x 96 wells
Format:	10 x 96-well
Assay Type:	Enzyme-Linked ImmunoSpot (ELISpot)
Signal:	Colorimetric
Specificity:	Recognizes natural murine IFN-g
Cross Reactivity:	No cross reactivity with other murine cytokines
Components:	ELISpot matched antibody pairs are extensively validated and include pre-titrated capture antibody and biotinylated detection antibody. Antibodies are supplied in quantities sufficient for 10 x 96 samples.
Background:	an immune molecule produced by T lymphocytes in response to mitogens or antigens [RGD, Feb 2006]. Publication Note: This RefSeq record includes a subset of the publications that are available for this gene. Please see the Gene record to access additional publications. ##Evidence-Data-START#Transcript exon combination :: AF010466.1 [ECO:0000332] RNAseq introns :: single sample supports all introns SRS369743 [ECO:0000348] ##Evidence-Data-END##
Gene Symbol:	Ifng
Gene ID:	25712



[View online »](#)