



Human Periostin/OSF-2 Immunoassay

Catalog Number: EA800097

For the quantitative determination of human Periostin/OSF-2 concentrations
in cell culture supernates, serum, and plasma.

For research use only. Not for use in diagnostic procedures.

MANUFACTURED AND DISTRIBUTED BY:

OriGene Technologies, Inc.

9620 Medical Center Drive Suite 200
Rockville, MD 20850, USA

TABLE OF CONTENTS

SECTION	PAGE
BACKGROUND.....	1
PRINCIPLE OF THE ASSAY.....	1
TECHNICAL HINTS AND LIMITATIONS.....	2
PRECAUTIONS.....	2
KIT COMPONENTS& STORAGE CONDITIONS.....	3
OTHER SUPPLIES REQUIRED BUT NOT SUPPLIED.....	4
SPECIMEN COLLECTION & STORAGE.....	4
REAGENTS PREPARATION.....	4
ASSAY PROCEDURE	6
CALCULATION OF RESULTS.....	6
PERFORMANCE CHARACTERISTICS.....	8
REFERENCES.....	10



BACKGROUND

Periostin (POSTN, PN, or osteoblast-specific factor OSF-2) is a protein that in humans is encoded by the POSTN gene. Periostin functions as a ligand for alpha-V/beta-3 and alpha-V/beta-5 integrins to support adhesion and migration of epithelial cells. Periostin plays a critical role in the development of cardiac valves and in degenerative valvular heart disease. While periostin usually is localized to the subendothelial layer in healthy heart valves, its levels are highly increased in infiltrated inflammatory cells and myofibroblasts in angiogenic areas in atherosclerotic and rheumatic valvular heart disease in humans. Periostin has also been shown to increase the secretion of matrix metalloproteinase from valvular interstitial cells, endothelial cells, and macrophages. It is thought that periostin plays a role in cardiac valve complex degeneration by inducing both angiogenesis and matrix metalloproteinase production.

PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. A monoclonal antibody specific for Periostin/OSF-2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any Periostin/OSF-2 present is captured by the coated antibody after incubation. Following extensive washing, a biotin-conjugate antibody specific for Periostin/OSF-2 is added to detect the captured Periostin/OSF-2 protein in sample. For signal development, horseradish peroxidase (HRP)-conjugated Streptavidin is added, followed by tetramethylbenzidine (TMB) reagent. Following a wash to remove any unbound combination, and enzyme conjugate is added to the wells. Solution containing sulfuric acid is used to stop color development and the color intensity which is proportional to the quantity of bound protein is measurable at 450nm.

Schematic diagram:





1. This ELISA should not be used beyond the expiration data on the kit label.
2. To avoid cross-contamination, use a fresh reagent reservoir and pipette tips for each step.
3. To ensure accurate results, some details, such as technique, plasticware and water sources should be emphasized.
4. A thorough and consistent wash technique is essential for proper assay performance.
5. A standard curve should be generated for each set of samples assayed.
6. It is recommended that all standards and samples be assayed in duplicate.
7. Avoid microbial contamination of reagents and buffers. Buffers containing protein should be made under aseptic conditions and be prepared fresh daily.
8. In order to ensure the accuracy of the results, the standard curve should be made every time.

PRECAUTIONS

The Stop Solution suggested for use with this kit is an acid solution. Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

KIT COMPONENTS & STORAGE CONDITIONS

PART	SIZE	STORAGE OF OPENED/ RECONSTITUTED MATERIAL
------	------	--



Microwell Plate - antibody coated 96-well Microplate (8 wells ×12 strips)	1 plate	Return unused wells to the foil pouch containing the desiccant pack. Reseal along entire edge of the zip-seal. May be stored for up to 1 month at 2 – 8°C**
Standard - lyophilized,4000 pg/ml upon reconstitution	2 vials	Aliquot and Store at -20°C** for six months
Concentrated Biotin-Conjugated antibody (100X) - 120 ul/vial	1 vial	Store at 2-8°C **for six months
Concentrated Streptavidin-HRP solution (100X) - 120 ul/vial	1 vial	Store at 2-8°C** for six months
Standard /sample Diluent - 16 ml/vial	1 bottle	Store at 2-8°C** for six months
Biotin-Conjugate antibody Diluent - 16 ml/vial	1 bottle	Store at 2-8°C** for six months
Streptavidin-HRP Diluent - 16 ml/vial	1 bottle	Store at 2-8°C** for six months
Wash Buffer Concentrate (20x) - 30 ml/vial	1 bottle	Store at 2-8°C** for six months
Substrate Solution - 12 ml/vial	1 bottle	Store at 2-8°C** for six months
Stop Solution - 12 ml/vial	1 bottle	Store at 2-8°C** for six months
Plate Cover Seals	4 pieces	

**Provided this is within the expiration date of the kit.

OTHER SUPPLIES REQUIRED BUT NOT SUPPLIED

1. Microplate reader capable of measuring absorbance at 450 nm.
2. Pipettes and pipette tips.
3. Deionized or distilled water.



4. Squirrt bottle, manifold dispenser, or automated microplate washer.
5. 500 mL graduated cylinder.
6. Human Periostin/OSF-2 controls (optional; available from Solarbio).

SPECIMEN COLLECTION & STORAGE

Cell Culture Supernates - Centrifuge cell culture media at 1000×g to remove debris. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Serum - Use a serum separator tube (SST) and allow samples to clot for 2 hours at room temperature or overnight at 2-8°C. Centrifuge at approximately for 15 minutes at 1000×g. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Plasma - Collect plasma using EDTA, heparin, or citrate as an anticoagulant. Centrifuge for 15 minutes at 1000×g within 30 minutes of collection. Assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

Note: The normal human serum or plasma samples are suggested to make a 1:2 dilution.

REAGENTS PREPARATION

1. **Temperature returning** - Bring all kit components and specimen to room temperature (20-25°C) before use.
2. **Wash Buffer** - Dilute 30mL of Wash Buffer Concentrate with 570mL of deionized or distilled water to prepare 200mL of Wash Buffer. If crystals have formed in the concentrate Wash Buffer, warm to room temperature and mix gently until the crystals have completely dissolved.
3. **Standard\Specimen** - Reconstitute the Standard with 1.0mL of deionized or distilled water. This reconstitution produces a stock solution of 4000 pg/mL. Allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to making dilutions. Pipette 500µL of Standard/Specimen Diluent into 2000pg/ml tube and the remaining tubes. Use the stock solution of 4000pg/mL to produce a 2-fold dilution series (below). Mix each tube thoroughly and change pipette tips between each transfer. The 4000 pg/mL standard serves as the high standard. The Standard/specimen Diluent serves as the zero standard (0 pg/mL).

***If you do not run out of re-melting standard, store it at -20°C. Diluted standard shall not be reused.**


4. **Working solution of Biotin-Conjugate anti-human Periostin/OSF-2 antibody:** Make a 1:100 dilution of the concentrated Biotin-Conjugate solution

with the Biotin-Conjugate antibody Diluent in a clean plastic tube.


***The working solution should be used within one day after dilution.**

- Working solution of Streptavidin-HRP:** Make a 1:100 dilution of the concentrated Streptavidin-HRP solution with the Streptavidin-HRP Diluent in a clean plastic tube.

***The working solution should be used within one day after dilution.**



ASSAY PROCEDURE



↓ Aspirate and wash 5 times

Add 100µl Substrate solution to each well, incubate 15 minutes, 37°C. Protect from light.

↓

Add 50µl Stop solution to each well. Read at 450nm within 30 minutes.

CALCULATION OF RESULTS

1. The standard curve is used to determine the amount of specimens.
2. First, average the duplicate readings for each standard, control, and sample. All O.D. values are subtracted by the mean value of blank control before result interpretation.
3. Construct a standard curve by reducing the data using computer software capable of generating a four parameter logistic (4-PL) curve-fit. As an alternative, construct a standard curve by plotting the mean absorbance for each standard on the y-axis against the concentration on the x-axis and draw a best fit curve through the points on the graph.
4. The data may be linearized by plotting the log of the Periostin/OSF-2 concentrations versus the log of the O.D. and the best fit line can be determined by regression analysis. This procedure will produce an adequate but less precise fit of the data. If samples have been diluted, the concentration read from the standard curve must be multiplied by the dilution factor.
5. This standard curve is provided for demonstration only. A standard curve should be generated for each set of samples assayed.

Typical data using the Periostin/OSF-2 ELISA

Standard (pg/ml)	OD.	OD.	Average	Corrected
0	0.032	0.035	0.034	-----
62.5	0.139	0.135	0.137	0.103
125	0.215	0.219	0.217	0.184
250	0.373	0.379	0.376	0.343
500	0.711	0.706	0.709	0.675



1000	1.339	1.337	1.338	1.305
2000	2.265	2.263	2.264	2.231
4000	3.303	3.305	3.304	3.271

Absorbance(450nm)

Human Periostin/OSF-2 Concentration(pg/ml)
Representative standard curve for Periostin/OSF-2 ELISA.

Performance Characteristics

SENSITIVITY: The minimum detectable dose was 62.5 pg/mL.

SPECIFICITY: This assay recognizes both natural and recombinant human Periostin/OSF-2. The factors listed below were prepared at 50ng/ml in Standard /sample Diluent and assayed for cross-reactivity and no significant cross-reactivity or interference was observed.

Factors assayed for cross-reactivity



Recombinant human	Recombinant mouse	Recombinant rat
β IG-H3	Periostin/OSF-2	
Integrin α V β 3		
Integrin α V β 5		

REPEATABILITY: The coefficient of variation of both intra-assay and inter-assay were less than 10%.

RECOVERY: The recovery of Periostin/OSF-2 spiked to three different levels in four samples throughout the range of the assay in various matrices was evaluated.

Recovery of Periostin/OSF-2 in two matrices

Sample Type	Average % of Expected Range (%)	Range (%)
Citrate plasma	97	92-102
Cell culture supernatants	104	96-112

LINEARITY: To assess the linearity of the assay, three samples were spiked with high concentrations of Periostin/OSF-2 in various matrices and diluted with the appropriate Sample Diluent to produce samples with values within the dynamic range of the assay.

Dilution ratio	Recovery (%)	Citrate plasma	Cell culture supernatants
1:2	Average% of Expected	92	101
	Range (%)	88-96	96-106
1:4	Average% of Expected	97	107

	Range (%)	93-101	102-112
--	-----------	--------	---------

REFERENCES

1. Takeshita S, et al. (1993) The Biochemical Journal. 294. 294 (1): 271–8.
2. Morra L, Moch H (Nov 2011). Virchows Archiv. 459 (5): 465–75.
3. Conway SJ, et al. (Apr 2014). Cellular and Molecular Life Sciences. 71 (7): 1279–88.
4. Braun N, et al. (2010) Peritoneal Dialysis International. 33 (5): 515–28.
5. Zhou W, et al. (Feb 2015). Nature Cell Biology. 17 (2): 170–82.
6. Tilman G, et al. (2007). Molecular Cancer. 6 (80): 80.

