

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for EUD4101**

## Sodium Iodide Symporter (SLC5A5) (C-term) Rabbit Polyclonal Antibody

## **Product data:**

Product Type:	Primary Antibodies
Applications:	IF, IHC
Recommended Dilution:	<ul> <li>Immunofluorescence.</li> <li>Immunohistochemistry on Frozen Sections.</li> <li>Working Dilution: 1/1200-1/2000 using FITC conjugated secondaries.</li> <li>Incubation Time: Overnight at 2-8°C.</li> <li>Positive Control: Frozen sections of Rat thyroid.</li> <li>Suitable for studying the sensitivity of thyroid tumors to iodine treatment.</li> </ul>
Reactivity:	Human, Porcine, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide from the C-terminus of Rat NIS
Specificity:	This antibody recognizes the Na+/l- Symporter. Absorption with 50-100 μg immunogen per ml diluted antiserum abolishes staining.
Formulation:	State: Serum State: Lyophilized Serum
<b>Reconstitution Method:</b>	Restore in 100 µl distilled water
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	solute carrier family 5 member 5
Database Link:	<u>Entrez Gene 6528 Human</u> <u>Q92911</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Sodium lodide Symporter (SLC5A5) (C-term) Rabbit Polyclonal Antibody – EUD4101
Background:	Human Sodium lodide Symporter (hNIS) is responsible for iodide concentrating ability within thyroid follicular cells. It is a membrane bound glycoprotein with 13 membrane spanning domains and 14 extramembranous domains. It may represent an autoantigen in thyroid. Located in the cell membrane of thyroid cells, NIS can allow Sodium and lodine flow across the membrane.
Synonyms:	Sodium-iodide cotransporter, Na+ / I- Symporter, NIS

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US