

Product datasheet for **GA100079**

Human PACAP receptor (ADCYAP1R1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ADCYAP1R1 CRISPRa kit - CRISPR gene activation of human ADCYAP receptor type I
Symbol:	ADCYAP1R1
UniProt ID:	P41586
Synonyms:	PAC1; PAC1R; PACAPR; PACAPRI
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)