

Product datasheet for **GA100087**

Human ADH4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ADH4 CRISPRa kit - CRISPR gene activation of human alcohol dehydrogenase 4 (class II), pi polypeptide
Symbol:	ADH4
UniProt ID:	P08319
Synonyms:	ADH-2; HEL-S-4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)