

Product datasheet for **GA100090**

Human ADH7 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | ADH7 CRISPRa kit - CRISPR gene activation of human alcohol dehydrogenase 7 (class IV), mu or sigma polypeptide |
| Symbol: | ADH7 |
| UniProt ID: | <u>P40394</u> |
| Synonyms: | ADH4 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)