

Product datasheet for **GA100096**

Human Adenosine A3 Receptor (ADORA3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ADORA3 CRISPRa kit - CRISPR gene activation of human adenosine A3 receptor
Symbol:	ADORA3
UniProt ID:	P33765
Synonyms:	A3AR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)