

## Product datasheet for **GA100108**

### Human beta 3 Adrenergic Receptor (ADRB3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	ADRB3 CRISPRa kit - CRISPR gene activation of human adrenoceptor beta 3
Symbol:	ADRB3
UniProt ID:	<a href="#">P13945</a>
Synonyms:	BETA3AR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)